


EXPERIENCE STRENGTH & STABILITY



Choice of 4-hole plate (30-38mm) and 2-hole plate (28-36mm)
ACCOMMODATES ANATOMICAL VERSATILITY

Large graft window
INCREASES INTRAOPERATIVE VISUALIZATION & POSITIONING

Intuitive single-step locking mechanism
FACILITATES SECURE PLACEMENT OF PLATE




Low profile, (3.5mm plate thickness) size matched plates
ALLOW FOR MULTI-LEVEL CONSTRUCTS


Multi-point fixation
PROVIDES OPTIMUM STABILITY

15° cephalad/caudal screw angulation
PRESENTS FLEXIBILITY IN PLATE SELECTION TO ASSIST IN FUSION

4-Hole Lateral Plate



2-Hole Lateral Plate



PRECISION SPINE®
Discover the Difference




SYSTEM TRAYS

P/N	Description
58-BK-1000	Tray for 4-Hole Plate
58-BK-1001	Tray for 2-Hole Plate




AccuFit® LATERAL PLATES & INSTRUMENTS

4-HOLE PLATES	
P/N	Description
58-LP-0030	AccuFit Lateral Plate - 30mm
58-LP-0032	AccuFit Lateral Plate - 32mm
58-LP-0034	AccuFit Lateral Plate - 34mm
58-LP-0036	AccuFit Lateral Plate - 36mm
58-LP-0038	AccuFit Lateral Plate - 38mm



2-HOLE PLATES	
P/N	Description
58-LP-0228	AccuFit Lateral Plate - 28mm
58-LP-0230	AccuFit Lateral Plate - 30mm
58-LP-0232	AccuFit Lateral Plate - 32mm
58-LP-0234	AccuFit Lateral Plate - 34mm
58-LP-0236	AccuFit Lateral Plate - 36mm



INSTRUMENTS	
P/N	Description
58-DA-0302	Angled Awl, Trocar Tip
58-DA-0301	Flexible Awl, Trocar Tip
58-DA-0201	Flexible Awl, Bevel Tip
58-DA-0101	Flexible Drill
58-LS-0100	Locking Screw
58-SD-0101	Flexible Screw Driver
58-RD-1375	7.5mm Up-Biting Rongeur
58-DA-0300	Awl, Trocar Tip
58-DA-0200	Awl, Bevel Tip
58-DA-0100	Drill
58-SD-0100	Screw Driver
58-HD-0100	Ratcheting AO Handle
58-SD-0400	Torque Limiting Handle
58-PH-0100	4-Hole Plate Holder
58-PH-0200	2-Hole Plate Holder
58-DG-0100	Variable Drill Guide
58-SD-0300	Lock/Unlock Driver
58-SD-0200	Lock Screw Driver

INSTRUMENTS BY REQUEST FOR 2-HOLE PLATES

P/N	Description
58-PH-0300	2-Hole 7° DTS/Guide Plate Holder
58-DA-0325	2-Hole DTS/Guide Drill – 25mm

AccuFit® LATERAL PLATE SCREWS

Ø 5.5mm VARIABLE SCREWS*	
P/N	Description
58-VS-5530	5.5 x 30mm Variable Screw
58-VS-5535	5.5 x 35mm Variable Screw
58-VS-5540	5.5 x 40mm Variable Screw
58-VS-5545	5.5 x 45mm Variable Screw
58-VS-5550	5.5 x 50mm Variable Screw
58-VS-5555	5.5 x 55mm Variable Screw
58-VS-5560	5.5 x 60mm Variable Screw

*25mm available by Special Request



Ø 6.5mm VARIABLE SCREWS*	
P/N	Description
58-VS-6530	6.5 x 30mm Variable Screw
58-VS-6535	6.5 x 35mm Variable Screw
58-VS-6540	6.5 x 40mm Variable Screw
58-VS-6545	6.5 x 45mm Variable Screw
58-VS-6550	6.5 x 50mm Variable Screw
58-VS-6555	6.5 x 55mm Variable Screw
58-VS-6560	6.5 x 60mm Variable Screw

*25mm available by Special Request

