


EXPERIENCE ZERO PROFILE



40° screw angle
RESISTS EXPULSION FORCES

Large bone graft window
ALLOWS FOR EXTENSIVE BONE GROWTH

Tantalum markers
PROVIDE EXCELLENT VISIBILITY AND POSITIONING



Four divergent screws (10°)
PROVIDE STRONG & STABLE FIXATION

Single-step locking screw
PREVENTS BONE SCREW BACKOUT

8° & 15° Lordotic options
ADDRESS DIFFERENT ANATOMICAL REQUIREMENTS

PRECISION SPINE®
Discover the Difference



SYSTEM TRAYS

P/N	Description
21-1040-CA-2	Vault Implant Tray
21-1041-CA-2	Vault Implant Tray



VAULT® PLATES & CAGES

8° - 32mm PLATES & CAGES

P/N	Description
10-P3208-11P	Vault Plate - 32mm x 11mm, 8°
10-P3208-13P	Vault Plate - 32mm x 13mm, 8°
10-P3208-15P	Vault Plate - 32mm x 15mm, 8°
10-P3208-17P	Vault Plate - 32mm x 17mm, 8°
10-P3208-19P	Vault Plate - 32mm x 19mm, 8°
10-C3208-11	Vault Cage - 32mm x 11mm, 8°
10-C3208-13	Vault Cage - 32mm x 13mm, 8°
10-C3208-15	Vault Cage - 32mm x 15mm, 8°
10-C3208-17	Vault Cage - 32mm x 17mm, 8°
10-C3208-19	Vault Cage - 32mm x 19mm, 8°

15° - 32mm PLATES & CAGES

P/N	Description
10-P3215-11P	Vault Plate - 32mm x 11mm, 15°
10-P3215-13P	Vault Plate - 32mm x 13mm, 15°
10-P3215-15P	Vault Plate - 32mm x 15mm, 15°
10-P3215-17P	Vault Plate - 32mm x 17mm, 15°
10-P3215-19P	Vault Plate - 32mm x 19mm, 15°
10-C3215-11	Vault Cage - 32mm x 11mm, 15°
10-C3215-13	Vault Cage - 32mm x 13mm, 15°
10-C3215-15	Vault Cage - 32mm x 15mm, 15°
10-C3215-17	Vault Cage - 32mm x 17mm, 15°
10-C3215-19	Vault Cage - 32mm x 19mm, 15°

8° - 39mm PLATES & CAGES

P/N	Description
10-P3908-11P	Vault Plate - 39mm x 11mm, 8°
10-P3908-13P	Vault Plate - 39mm x 13mm, 8°
10-P3908-15P	Vault Plate - 39mm x 15mm, 8°
10-P3908-17P	Vault Plate - 39mm x 17mm, 8°
10-P3908-19P	Vault Plate - 39mm x 19mm, 8°
10-C3908-11	Vault Cage - 39mm x 11mm, 8°
10-C3908-13	Vault Cage - 39mm x 13mm, 8°
10-C3908-15	Vault Cage - 39mm x 15mm, 8°
10-C3908-17	Vault Cage - 39mm x 17mm, 8°
10-C3908-19	Vault Cage - 39mm x 19mm, 8°

15° - 39mm PLATES & CAGES

P/N	Description
10-P3915-11P	Vault Plate - 39mm x 11mm, 15°
10-P3915-13P	Vault Plate - 39mm x 13mm, 15°
10-P3915-15P	Vault Plate - 39mm x 15mm, 15°
10-P3915-17P	Vault Plate - 39mm x 17mm, 15°
10-P3915-19P	Vault Plate - 39mm x 19mm, 15°
10-C3915-11	Vault Cage - 39mm x 11mm, 15°
10-C3915-13	Vault Cage - 39mm x 13mm, 15°
10-C3915-15	Vault Cage - 39mm x 15mm, 15°
10-C3915-17	Vault Cage - 39mm x 17mm, 15°
10-C3915-19	Vault Cage - 39mm x 19mm, 15°

VAULT® SCREWS & LOCKING CAPS


LOCKING CAPS

P/N	Description
10-1002-11	Vault Locking Cap - 11mm
10-1002-13	Vault Locking Cap - 13mm
10-1002-15	Vault Locking Cap - 15mm
10-1002-17	Vault Locking Cap - 17mm
10-1002-19	Vault Locking Cap - 19mm



Ø 4.0mm SELF-DRILLING VARIABLE SCREWS

P/N	Description
10-4025	Vault SDV Screw - 4.0mm x 25mm
10-4030	Vault SDV Screw - 4.0mm x 30mm
10-4035	Vault SDV Screw - 4.0mm x 35mm



SET SCREW

P/N	Description
10-1001	Vault Set Screw

