

EXPERIENCE ZERO PROFILE



42° screw angulation
MAXIMIZES BONE PURCHASE & STABILITY


Large graft window
OPTIMIZES BIOLOGICAL
INGROWTH TO
FACILITATE FUSION

One-step locking mechanism
PREVENTS BONE SCREW BACKOUT

Ti plate & integrated PEEK cage
PROVIDE ANTERIOR COLUMN
SUPPORT & STABILITY

Three footprints
ACCOMMODATE VARYING PATIENT ANATOMY

0° & 7° Lordotic options
ASSIST IN REPRODUCING
PATIENT'S SAGITTAL PROFILE



PRECISION SPINE®
Discover the Difference



SYSTEM TRAYS

P/N	Description
37-BK-01XX	Vault C Standard Implants & Sizers
37-BK-021X	Vault C Rasps (By request)
37-BK-0201	Vault C Instruments



18x15mm footprint available by request



CAGES & PLATES

14x12mm 0° CAGES

P/N	Description
37-CP-4206	Vault C Cage 14x12 - 0° x 6mm
37-CP-4207	Vault C Cage 14x12 - 0° x 7mm
37-CP-4208	Vault C Cage 14x12 - 0° x 8mm
37-CP-4209	Vault C Cage 14x12 - 0° x 9mm
37-CP-4210	Vault C Cage 14x12 - 0° x 10mm
37-CP-4211	Vault C Cage 14x12 - 0° x 11mm
37-CP-4212	Vault C Cage 14x12 - 0° x 12mm

14x12mm 7° CAGES

P/N	Description
37-CL-4206	Vault C Cage 14x12 - 7° x 6mm
37-CL-4207	Vault C Cage 14x12 - 7° x 7mm
37-CL-4208	Vault C Cage 14x12 - 7° x 8mm
37-CL-4209	Vault C Cage 14x12 - 7° x 9mm
37-CL-4210	Vault C Cage 14x12 - 7° x 10mm
37-CL-4211	Vault C Cage 14x12 - 7° x 11mm
37-CL-4212	Vault C Cage 14x12 - 7° x 12mm

14x12mm PLATES

P/N	Description
37-PA-4206	Vault C Plate 14x12 - 6mm
37-PA-4207	Vault C Plate 14x12 - 7mm
37-PA-4208	Vault C Plate 14x12 - 8mm
37-PA-4209	Vault C Plate 14x12 - 9mm
37-PA-4210	Vault C Plate 14x12 - 10mm
37-PA-4211	Vault C Plate 14x12 - 11mm
37-PA-4212	Vault C Plate 14x12 - 12mm

16x14mm 0° CAGES

P/N	Description
37-CP-6406	Vault C Cage 16x14 - 0° x 6mm
37-CP-6407	Vault C Cage 16x14 - 0° x 7mm
37-CP-6408	Vault C Cage 16x14 - 0° x 8mm
37-CP-6409	Vault C Cage 16x14 - 0° x 9mm
37-CP-6410	Vault C Cage 16x14 - 0° x 10mm
37-CP-6411	Vault C Cage 16x14 - 0° x 11mm
37-CP-6412	Vault C Cage 16x14 - 0° x 12mm

16x14mm 7° CAGES

P/N	Description
37-CL-6406	Vault C Cage 16x14 - 7° x 6mm
37-CL-6407	Vault C Cage 16x14 - 7° x 7mm
37-CL-6408	Vault C Cage 16x14 - 7° x 8mm
37-CL-6409	Vault C Cage 16x14 - 7° x 9mm
37-CL-6410	Vault C Cage 16x14 - 7° x 10mm
37-CL-6411	Vault C Cage 16x14 - 7° x 11mm
37-CL-6412	Vault C Cage 16x14 - 7° x 12mm

16x14mm PLATES

P/N	Description
37-PA-6406	Vault C Plate 16x14 - 6mm
37-PA-6407	Vault C Plate 16x14 - 7mm
37-PA-6408	Vault C Plate 16x14 - 8mm
37-PA-6409	Vault C Plate 16x14 - 9mm
37-PA-6410	Vault C Plate 16x14 - 10mm
37-PA-6411	Vault C Plate 16x14 - 11mm
37-PA-6412	Vault C Plate 16x14 - 12mm

18x15mm 0° CAGES

P/N	Description
37-CP-8506	Vault C Cage 18x15 - 0° x 6mm
37-CP-8507	Vault C Cage 18x15 - 0° x 7mm
37-CP-8508	Vault C Cage 18x15 - 0° x 8mm
37-CP-8509	Vault C Cage 18x15 - 0° x 9mm
37-CP-8510	Vault C Cage 18x15 - 0° x 10mm
37-CP-8511	Vault C Cage 18x15 - 0° x 11mm
37-CP-8512	Vault C Cage 18x15 - 0° x 12mm

18x15mm 7° CAGES

P/N	Description
37-CL-8506	Vault C Cage 18x15 - 7° x 6mm
37-CL-8507	Vault C Cage 18x15 - 7° x 7mm
37-CL-8508	Vault C Cage 18x15 - 7° x 8mm
37-CL-8509	Vault C Cage 18x15 - 7° x 9mm
37-CL-8510	Vault C Cage 18x15 - 7° x 10mm
37-CL-8511	Vault C Cage 18x15 - 7° x 11mm
37-CL-8512	Vault C Cage 18x15 - 7° x 12mm

18x15mm PLATES

P/N	Description
37-PA-8506	Vault C Plate 18x15 - 6mm
37-PA-8507	Vault C Plate 18x15 - 7mm
37-PA-8508	Vault C Plate 18x15 - 8mm
37-PA-8509	Vault C Plate 18x15 - 9mm
37-PA-8510	Vault C Plate 18x15 - 10mm
37-PA-8511	Vault C Plate 18x15 - 11mm
37-PA-8512	Vault C Plate 18x15 - 12mm

SCREWS

Ø 3.5mm SELF-DRILLING SCREWS

P/N	Description
37-SD-3512	SD Screw - 3.5 x 12mm
37-SD-3514	SD Screw - 3.5 x 14mm
37-SD-3515	SD Screw - 3.5 x 15mm
37-SD-3517	SD Screw - 3.5 x 17mm

Ø 3.5mm BLUNT TIP, SELF-TAPPING SCREWS

P/N	Description
37-SB-3512	ST Screw - 3.5 x 12mm
37-SB-3514	ST Screw - 3.5 x 14mm
37-SB-3515	ST Screw - 3.5 x 15mm
37-SB-3517	ST Screw - 3.5 x 17mm

Ø 4.0mm SELF-DRILLING SCREWS

P/N	Description
37-SD-4012	SD Screw - 4.0 x 12mm
37-SD-4014	SD Screw - 4.0 x 14mm
37-SD-4015	SD Screw - 4.0 x 15mm
37-SD-4017	SD Screw - 4.0 x 17mm

Ø 4.0mm BLUNT TIP, SELF-TAPPING SCREWS

P/N	Description
37-SB-4012	ST Screw - 4.0 x 12mm
37-SB-4014	ST Screw - 4.0 x 14mm
37-SB-4015	ST Screw - 4.0 x 15mm
37-SB-4017	ST Screw - 4.0 x 17mm

18x15mm footprint available by request

PRECISION SPINE®
www.precisionspineinc.com

